
THE ANAPEN® TRAINING DEVICE

The Anapen® Training Device (Trainer) is a dummy device which allows the user to become familiar with the device and how it works. **The Trainer does not have a needle and does not contain adrenaline.** Anapen® Trainers can be used repeatedly but will tire with time and excessive use. Follow the instructions for firing and resetting to maximise the life of the device.



FIRING THE ANAPEN® TRAINER



Remove the black needle shield cap A by pulling hard, in the direction of the arrow. This also removes a grey protective needle shield ©.



Remove the grey safety cap (B) from the red firing button (C), by pulling as indicated by the arrow.



Hold the open end of the device (E) against the outer part of your thigh.



Once the open end of the device (s) is against the outer part of your thigh, press the red firing button (s) so that it clicks. Keep holding the Anapen® Auto-Injector against your thigh for 3 seconds. Slowly remove Anapen® from your thigh.




The injection indicator
will have turned red. This shows that the injection is complete (If the injection indicator is not red, you must repeat the injection with a new Anapen®.)



After the injection, click the wide end of the black needle shield cap (A) back on the open end of the (E) device (as indicated by the arrow).

HOW TO RESET THE ANAPEN® TRAINER



Replace the grey protective needle shield (G) on dummy syringe (F).

Before moving to the next step, place the grey protective needle shield back onto the end of the plastic dummy syringe.



Fit the open end of the device $\[\]$ into the thin end of the black needle shield cap. $\[\]$



Push the black cap (A) until it clicks with the help of a hard surface.



Replace the grey safety cap (B) on the red firing button.

Stand the training pen upright with the black boot remover against a hard surface, and reattach the grey safety cap.

Anapen® TRAINER is now ready for practice use again.

Distributed by Allergy Concepts Pty Ltd Level 1, 83 York Street Sydney NSW 2000 Email: info@anapen.com.au www.anapen.com.au NOV 2022

