

URINSAFE™

Reagent Strips for Urinalysis

| pH | Meaning |
|------------|--------------------------------|
| 5.0 to 7.0 | Normal |
| 7.0 to 8.0 | Normal - add .005 to SG result |
| 8.0 | Caution* |

| Specific Gravity | Hydration Status |
|------------------|-----------------------|
| 1.000 to 1.005 | Excellent ** |
| 1.005 to 1.010 | Very Good |
| 1.010 to 1.015 | Good |
| 1.015 to 1.020 | Fair |
| 1.020 to 1.025 | Caution - Borderline |
| >1.025 | Poor - Unsatisfactory |

100 Strips


TESTING AND READING TIME

Rev. 04/2019

| | | | | | | | |
|---------------------------------|-------|-------|-------|-------|-------|-------|-------|
| pH 60s | 5.0 | 6.0 | 6.5 | 7.0 | 7.5 | 8.0 | 8.5 |
| Specific Gravity 45s | 1.000 | 1.005 | 1.010 | 1.015 | 1.020 | 1.025 | 1.030 |



IVD



* A pH of 8.0 **may** indicate a urinary tract infection or an old sample. Retest with a fresh sample. If the pH is still 8.0 seek medical advice.

** With persistent low readings, **approaching an SG of 1.000**, seek medical advice as it is necessary to rule out over hydration or dilutional hyponatremia.

LOT



Distributed by
LFA First Response
 Australian helpline 1800 681 544
www.lfafirstresponse.com.au