

### X SERIES® ADVANCED

The new best in class monitor/defibrillator from ZOLL



ZOLL's new X Series® Advanced monitor/defibrillator helps providers manage patients more effectively than ever before. X Series Advanced is the first all-in-one device of its kind to provide the clinical data EMS teams need at every step to improve outcomes.



### Real BVM Help™ and TBI Dashboard

With groundbreaking new technology available only through ZOLL, X Series Advance finally addresses some of the key challenges providers face in both the pre-hospital and hospital environments today.



## Real BVM Help – Real-time feedback for ventilation

With Real BVM Help on the X Series Advanced, providers can view delivered tidal volume and rate directly on the monitor. This real-time feedback guides users to deliver high-quality ventilation in a manner that has never been possible – until now.

### Real BVM Help on the ZOLL X Series Advanced gives providers:

- Target and delivered ventilation volume
- Target and delivered ventilation rate
- A ventilation quality indicator and countdown timer



# TBI Dashboard: The information you need when it matters most

Traumatic Brain Injury (TBI) patients present a unique set of challenges. Early treatment is critical, as additional secondary injury often occurs after the initial event.<sup>1</sup>


ZOLL's TBI Dashboard on the X Series Advanced monitor/defibrillator provides trending information on the parameters most critical to a TBI patient, enabling clinicians to quickly detect potential patient deterioration.

<sup>1</sup>Spaite DW, et al. Acad Emerg Med. 2014;21:7:818-830



#### Meet your telemedicine needs with Remote View on the X Series Advanced

With Remote View, clinicians can view data from the X Series Advanced from a remote location. This telemedicine solution from ZOLL enables effective decision support to allow EMS providers meet their telehealth needs



#### Real BVM Help integrates with existing ZOLL CPR technologies

For cardiac arrest patients, the X Series Advanced monitor/defibrillator offers both Real CPR Help® and Real

BVM Help, giving providers simultaneous real-time feedback on both CPR and ventilation quality. When used with adult ZOLL CPR electrodes, the X Series Advanced device displays the Real BVM Help Dashboard next to the Real CPR Help Dashboard.





## Improve your system's performance with RescueNet CaseReview

RescueNet® CaseReview gives providers the ability to view data from the X Series Advanced post-case for effective QA/QI. Data available in CaseReview includes ventilation performance with Real BVM Help, CPR performance with Real CPR Help®, and automatic recognition and integration of performance with ZOLL's automated CPR device, the AutoPulse®.

### X Series Advanced: The new best in class monitor/defibrillator from ZOLL.

For more information contact **LFA First Response**:

www.lfafirstresponse.com.au | 1800 681 544

