



LFA135



PHARMA 1000

UNDER COUNTER, OR COUNTER TOP VACCINE REFRIGERATOR

SPECIFICATIONS



TEMPERATURE RANGE

+3°C to +7°C
@ 10°C to 32°C room temperature
& 60% humidity



CAPACITY

Effective capacity: 145 Litres



DIMENSIONS

Width: 595mm
Depth: 640mm
Height: 830mm

FEATURES



FEATURES

- Purpose built for vaccine storage
- Exceptional performance on all shelves
- Full use of internal space (except 4cm below the fan)
- QCPP approved and exceeds requirements
- No frozen back wall
- Simple to use and no setup required
- Very easy to use hi/lo/reset temperature display
- Can be built in – only a 3cm top gap is required
- Long lasting LED lighting
- Available with right or left hinged door
- Designed and manufactured in Denmark
- 4 adjustable shelves plus raised base shelf
- Lockable and self-closing glass door
- Energy consumption – only 0.7 kWh/24h
- Extremely quiet – only 35dB
- Audible and visual alarm
- Voltage free external alarm
- Data logger access valve
- Power failure alert for Building Maintenance System
- Built in door alarm