stryker



Stair-PRO®

stair chair

As a medic you are relentless. And in your relentless pursuit for optimal patient care, you are faced with unique challenges patient handling in stairways is a challenge you may face on any given call. The strain on your back from lifting and lowering patients down stairs can add up over time. With 1 in 4 EMS workers suffering a career-ending back injury within the first four years of employment, it is clear you need an option to reduce this risk¹.

Each patient handled in a stairway can pose a risk for caregiver injury. Stair-PRO puts the power to help reduce that risk back in your hands. Stryker's Stair-PRO stair chair – it's a powerful thing.



Compact folded depth

Easy-to-fold functionality allows the chair to fold to a compact size for convenient storage in existing vehicle compartments.

Optional accessories:







Foot Support



IV Clip



Vinyl head support

Key features:



Specifications¹

• Models	6252	
Height	37.5 in	(95.25 cm)
• Width	20.5 in	(52 cm)
• Depth ²	28 in	(71 cm)
Folded Depth	8 in	(20 cm)
• Weight	32.8 lb	(14.9 kg)
Maximum Load ³	500 lb	(228 kg)

¹ Dimensions are measured from the outermost edges of the main frame. Speci cations are rounded to the nearest whole number. Conversions are calculated before rounding.

In-service video included with every order.

Stryker reserves the right to change speci cations without notice.

References

1. Sanders, Mick J. (2011) Mosby's Paramedic Textbook (4th ed., p. 36)

Certifications



Waranty:

- One-year parts, labor and travel or two-year parts only
- Lifetime on all welds*

Extended waranties available:

• 7-year service life

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Flex Financial, ProCare, Stair-PRO, Stryker. All other trademarks are trademarks of their respective owners or holder.

The yellow and black color sceme is a registered trademark of Stryker Corporation

Technical Support

Stryker Technical Support comprises a team of professionals available to help with your Stair-PRO needs. Contact via phone at 1 800 STRYKER or email at ${\tt Technical Support@Stryker.com}.$

ProCare Service from Stryker

With maintenance, service and work-flow solutions that proactively maintain your equipment, we help you get every last penny out of your investment. Two-hour priority response callback time and 24- to 72-hour product turnaround help you focus more on your patients and less on equipment upkeep.

Flex Financial Program

Our financial programs provide a range of smart alternatives designed to fit your organization's needs. We offer flexibility beyond a cash purchase with payment structures customizable to meet budgetary needs and help to build long-term financial stability. Specialized finance solutions for municipalities and other government agency businesses are also available.

To learn more about the wide range of financial services, call 1 888 308 3146.

3800 E. Centre Avenue Portage, MI 49002 USA t: 269 329 2100

toll free: 800 327 0770

ems.stryker.com

 $^{^{2}}$ Depth dimensions are measured with extendable handles retracted.

³ Maximum load capacity is total weight distributed in accordance to basic human anatomy. EMTs must consider the weight of the patient, equipment and accessories when determining the total weight on the product.

⁴Patient restraint set includes two cross-chest restraints and one ankle restraint. Cross-chest restraints can also be used in a chest and lap restraint con guration.