



VivaGuard™ FT-100A
Infrared Forehead Thermometer



CE 0197



Quick start guide



Make sure the forehead not covered with hair, sweat or dirt



Aim at forehead with distance within 2 cm



Press the **SCAN** button, then read the result



1
Sec

1 Second Measurement

°C / °F

°C / °F Switchable

20
MEM

20 Memory



Red backlight fever alert

Specification

Analyzer	Specification
Range	32.0 - 43.0°C (90.0 - 109.4°F)
Accuracy	± 0.2°C (± 0.4°F)
°F/°C switchable	Yes
Jumbo LCD	Yes
Waterproof	No
Fever Alarm	Yes
Battery Life	> 3,000 tests
Battery replaceable	Yes
Auto-Shut off	Yes
Backlight	Yes
Beep function	Yes
Measurement modes	Body/Object mode switchable
Sensor type	Infrared sensor
Thermometer dimension	162 x 50 x 35 mm
Weight	About 98g (with batteries)
Bluetooth	Optional
Resolution	0.1°C /°F
Warranty	2 years