



PHARMA 3000

FLOOR STANDING VACCINE REFRIGERATOR



TEMPERATURE RANGE

+3°C to +7°C @ 10°C to 32°C room temperature & 60% humidity



CAPACITY

Effective capacity: 350 Litres



DIMENSIONS

Width: 595mm

Depth: 640mm Height: 1900mm



FEATURES

- Purpose built for vaccine storage
- Exceptional performance on all shelves
- Full use of internal space (except 4cm below the fan)
- QCPP approved and exceeds requirements
- No frozen back wall and no cold spots
- Simple to use and no setup required
- Very easy to use hi/lo/reset temperature display
- Can be built in only a 3cm top gap is required
- Long lasting LED lighting (glass door model)
- 9 adjustable shelves plus raised base shelf

- Designed and manufactured in Denmark
- Lockable and self-closing, glass or solid door
- Available with right or left hinged door
- Energy consumption only 1.1 kWh/24h
- Extremely quiet only 36dB
- Audible and visual alarm
- Voltage free external alarm contacts for BMS
- Data logger access valve
- Environmentally responsible refrigerant
- Built in door alarm

