

BRAVDEN

Manikin Features



Oral and nasal passages which allow realistic nose pinch



Head tilt and chin lift for opening airway



Audible feedback reinforces correct compression depth (Clicker ON/OFF)



Highly realistic chest movement with correct ventilations



Anatomically correct. Harmless, latex-free skin. Easy AED training pad removal. RoHS & REACH compliant




Air vent reduces the risk of cross contamination. Air from Brayden's disposable lungs exits through the back of his head. This means students never come in to contact with exhaled air

DID YOU KNOW...

Brayden has been tested to endure over 1,000,000 compressions



Brayden Kit Contents



DID YOU KNOW...

Brayden uses state-of-the-art technology with fewer moving parts and less wires resulting in a much longer life

