

Devac

For ultra clean suction lines

- Easy to use
- Two-part system
- +
- Non-foaming
- +
- Removes biofilm



- 💧 Foam-free liquid concentrate for professional cleaning, deodorisation and protecting all suction systems and spittoons.
- 💧 Removes conventional biological soils and bacterial biofilms.
- 💧 Regular use of Devac DAILY & WEEKLY will increase longevity and life span of suction systems.



Try DEVAC now for FREE when you sign up for our 2 week trial.



DEVAC 2 WK TRIAL SCHEDULE	MON	TUES	WEDS	THURS	FRI
Week 1	Use 25mL Devac Daily	Use 25mL Devac Daily	Use 25mL Devac Daily	Use 25mL Devac Daily	Use 25mL Devac Weekly
Week 2	Use 25mL Devac Daily	Use 25mL Devac Daily	Use 25mL Devac Daily	Use 25mL Devac Daily	Use 25mL Devac Weekly

Devac

DAILY

Devac DAILY is highly effective in removing debris and bio-burden from evacuation lines and pumps. It provides excellent maintenance and ensures optimal suction on all dental aspiration systems. The specially formulated alkaline based formulation is non-foaming. This product is used as a 2.5% solution [25ml of product to 1 litre of water].

	PACKSIZE	CODE
Devac DAILY	6 x 1L	040232

Devac

WEEKLY

Devac WEEKLY is designed to provide heavy duty removal of debris and bio-burden from evacuation lines and pumps. Used regularly it provides excellent maintenance and ensures optimal suction on all dental aspiration systems. The specially formulated acid-based formulation is non-foaming. This product is used as a 2.5% solution [25ml of product to 1 litre of water].

	PACKSIZE	CODE
Devac WEEKLY	6 x 1L	040235

To claim your **FREE Trial** email marketing@whiteley.com.au today!