

HOW TO CLEAN SURFACES USING SPEEDY CLEAN BIODEGRADABLE WIPES



Hospital Grade Neutral Detergent Wipes



Made from 100% Plant Based Fibre



Plastic Free Wipes



100% Dye Free



Wipe fibres will biodegrade under certified commercial composting conditions

Additional information listed in the SDS – whiteley.com.au



ARTG No: 125529



01

Personal Protective Equipment

Ensure appropriate PPE is worn prior to commencing.



02

Clean and remove soiling

Remove gross soiling from the surface and clean using a Speedy Clean Biodegradable Wipe in one direction only across the surface to prevent contamination. Dispose of the soiled wipe to an appropriate receptacle.



03

Repeat as required

Repeat as required until surface is visibly clean. Allow surface to dry.